ASEPTIQUIK® S SERIES CONNECTORS



AseptiQuik® S Small Format Connectors provide quick and easy sterile connections for small-flow applications, even in non-sterile environments. The "FLIP-CLICK-PULL" design of AseptiQuik S enables users to quickly transfer small volumes of media easily with less risk of operator error. The connector's genderless and robust design provides reliable performance without the need for clamps, fixtures or tube welders. Biopharmaceutical manufacturers can now make 1/8", 1/4" and 3/8" hose barb and 1/4" and 3/4" sanitary sterile connections with the quality and market availability they expect from the leader in single-use connection technology.

FEATURES BENEFITS Genderless design Eases integration of single-use systems with one part number for both halves Robust construction Repeatable and reliable performance with no additional hardware required FLIP-CLICK-PULL Innovative three-step connection process reduces risk of operator error Integrated pull tab covers Pull tabs act as protective cover reducing part complexity and ensure simultaneous removal of both membranes

Audible confirmation of assembly

throughout upstream processing

Minimizes transitional flow disruptions

Specifications • • •







PRESSURE:

Up to 60 psi, 4.1 bar

TEMPERATURE:

39°F to 104°F (4°C to 40°C)

STERILIZATION:

Gamma: Up to 50kGy irradiation Autoclave High Temp (HT) Version: Up to 266°F (130°C) for 60 minutes

TERMINATION SIZES:

1/8", 1/4" and 3/8" ID hose barb (3.2mm, 6.4mm and 9.5mm), 1/4" and 3/4" sanitary and MPC insert

MATERIALS:

Main Components:

Polycarbonate (white), USP Class VI, ADCF

Pull Tabs/Caps:


Polycarbonate (blue, standard version), USP Class VI, ADCF Polycarbonate (white, HT version), USP Class VI, ADCF

Silicone (clear), platinum-cured, USP Class VI, ADCF

Membrane:

Polyethylene (standard version), USP Class VI, ADCF Hydrophobic polyethersulfone (HT version), USP Class VI, ADCF, PTFE strip sticker

This graph is intended to give you a general idea of the performance capabilities of each product line. The shaded area of the graph represents the operating range of the product family, i.e. upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area. AseptiQuik S Flow



Visit us at cpcworldwide.com or call 800-519-7633 to get a free catalog or to find your local distributor.

CPC Click

1/4 sanitary termination


with straight bore


ASEPTIQUIK® S SERIES DIMENSIONS

A = Height/Diameter

B = Total Length

Hose Barb Length





Coupling Bodies

POLYCARBONATE with Blue Pull Tabs



TERMINATION

1/8" HOSE BARB

1/4" HOSE BARB

1/4" SANITARY

3/8" HOSE BARB

3/4" SANITARY

MPC INSERT

PART NO.
AQS17002
AQS17004
AQS33004
AQS17006
AQS33012
AQS17MPC

2.25 2.25 2.25 2.25 2.25 2.25 2.25

1.30 1.45 1.50 1.45 1.60 1.49 0.50 0.65 0.70 0.65 0.80 0.69

POLYCARBONATE HT with White Pull Tabs



TERMINATION

1/8" HOSE BARB

1/4" HOSE BARB

1/4" SANITARY

3/8" HOSE BARB

3/4" SANITARY

MPC INSERT

PART NO.
AQS17002HT
AQS17004HT
AQS33004HT
AQS17006HT
AQS33012HT
AQS17MPCHT

B 1.30 1.45 1.50 1.45 1.60 1.49 0.50 0.65 0.70 0.65 0.80 0.69

Flip-Click-Pull Assembly Procedure



Unsnap and flip down the protective pull tab covers on each connector half.



2



Align the connector halves with the pull tabs hanging down. Slide the two halves together. Independently squeeze each side of the connector until you hear an audible "CPC Click."



To complete the connection, simply snap the pull tabs together by pushing on the CPC logos and pull the membranes from the connector.



COLDER PRODUCTS COMPANY

PHONE +1 651-645-0091 FAX +1 651-645-5404 TOLL FREE 800-519-7633 COLDER PRODUCTS COMPANY GMBH
Germany

PHONE +49-6026-9973-0 FAX +49-6026-9973-173 COLDER PRODUCTS COMPANY LIMITED Hong Kong S.A.R. of China

PHONE +852-2987-5272 FAX +852-2987-2509

Confidence at every point of connection.

cpcworldwide.com/ContactUs

WARRANTY: All sales are subject to Colder Products Company's limited express warranty set forth in the CPC catalog. Contact your local distributor or CPC Customer Service for warranty provisions. Warning: Due to the wide variety of possible fluid media and operating conditions, unintended consequences may result from the use of this product, all of which are beyond the control of CPC. It is the user's responsibility to carefully determine and test for compatibility for use with their application. All such risks shall be assumed by the buyer.