

## **DISPOSABLE SAMPLE VALVE**

SINGLE USE / LIMITED USE

Fluid Line Technology is proud to introduce our new *Disposable Sampling Valves*.

These disposable valves come in two different styles: **Single Use** and **Limited Use**. Single Use is designed to last 3 cleaning cycles while the Limited Use is designed for up to 12 cleaning cycles\*. This allows our customers to fine-tune their disposable process needs to minimize cost while maximizing throughput.

Each design is made from only certified USP CLASS IV materials with a steam sterilizing capability of 121°C for 30 minutes.

Each *Disposable Sampling Valve* comes in a sealed pouch (ready to sanitize) or a gamma irradiated pouch (ready to use). Contact us for more information about sizes and configuration options.

\* Refers to CIP or autoclave



	Features and Benefits	Single-Use	Limited-Use
>	4 cleaning cycles	•	
>	15 cleaning cycles		•
>	Fully autoclavable	•	•
>	USP Class VI materials	•	•
>	Ready-to-sanitize option	•	•
>	Gamma irradiation option	•	•
>	Adjustable flow control	•	•
>	Easy-to-use twist handle	•	•
>	TFM™ tip (seat)	•	•
>	Self-aligning & non-rotating stem	•	•
>	Zero dead leg for all sizes and connections**	•	•
>	No tools needed to operate	•	•

<sup>\*\*</sup>The internal dead volume is minimized to its fluid surface tension. This is not always true for our competitors TFM™ is registered by Dyneon, Inc.